

Oryca

Medical Supplies LLC

“Committed to Health, Driven by Quality”

Drug of Abuse Testing Solutions



(we recommend reading the results at 5-10 minutes)


Saliva Test Cube



Key Features

Format/Shape	cup
Time to result	at 5 minutes
Sample	saliva
Accuracy	over 99%
Certifications	CE, ISO13485
Customization	upto 28 customizable parameters


Urine Test Cassette



Key Features

Format/Shape	cassette
Time to result	at 5 minutes
Sample	urine
Accuracy	over 99%
Certifications	CLIA/510K, CE , ISO13485, FDA
Customization	upto 15 customizable parameters

Urine Dip Card



Key Features

Format/Shape	card
Time to result	at 5 minutes
Sample	urine
Accuracy	over 99%
Certifications	CLIA/510K, CE , ISO13485, FDA
Customization	upto 18 customizable parameters

Urine Round Cup



Key Features

Format/Shape	cup
Time to result	at 5 minutes
Sample	urine
Accuracy	over 99%
Certifications	CLIA/510K, CE , ISO13485, FDA
Customization	upto 24 customizable parameters with flexible cut-off levels

Round Cup Reader



Key Features

Time to result	cup scan in less than 40 seconds
Input power	15W
Size	225mm x 100mm x 160mm
Measurement range	4.0%~15.0%(20~140mmol/mol)
Weight	650g
Data export	through USB thumb drive or PC-link
Data save	saved as .csv & .pdf format
Additional Features	includes a barcode scanner & printer
Certifications	CE(IVDR)

Drug Testing Parameters and Cut-Off Levels

Test Name	Calibrator	Urine (ng/ml)	Oral Fluid (ng/ml)
Amphetamine/AMP	D-Amphetamine	100 / 300 / 500 / 1000	50
Acetaminophen/APAP	Acetaminophen	2000 / 5000	
Barbiturates/BAR	Secobarbital	200 / 300	50
Benzodiazepines/BZO	Oxazepam	5/10/ 50 / 100 / 150 / 200/ 300	10 / 50
Buprenorphine/BUP	Buprenorphine	5 / 10	10
Caffeine/CAF	Caffeine	8000	
Carfentanyl/CFYL	Carfentanyl	500	
Klonopin/CLZP	Klonopin	300	
Cocaine/COC	Benzoyllecgonine	10 / 50 / 100 / 150 / 200 / 300	20 / 50 / 150
Cotinine/COT	Cotinine	200	
EDDP/EDDP	EDDP	100 / 300	
Ethyl Glucuronide/ETG	Ethyl Glucuronide	300 / 500 / 1000	50
Fentanyl/FEN	Fentanyl	1 / 2 / 10 / 20 / 50 / 100 / 200 / 300	50
Gabapentin/GBPT	Gabapentin	1000 / 2000	
Hydrocodone/HDCD	Hydrocodone	100	
K2/K2	JWH-018/JWH-073	20 / 25 / 50	50
AB-PINACA/K3	AB-PINACA	25 / 300	100
K4/UR-144	UR-1445-Pentanoicacid Metabolite	10	
Ketamine/KET	Ketamine	100 / 300 / 500 / 1000	100
Kratom/KRA	Mitragynine	25 / 100 / 250 / 300 / 500	
Lysergic Acid Diethylamide/LSD	LSD	3 / 10	
Methcathinone/MCAT	S(-)-Methcathinone Hcl	1000	
Methadone/MTD	Methadone	50 / 150 / 200 / 300 / 500	35
Methamphetamine/mAMP/MET	D-Methamphetamine	100 / 300 / 500 / 1000	50
Methylenedioxyamphetamine/MDMA	MDMA	100 / 300 / 500	50 / 60
Methylenedioxypropylvalerone/MDPV	Methylenedioxypropylvalerone	300 / 500	
Mephedrone/MEP	Mephedrone	500	
Methaqualone/MQL	Methaqualone	300	50
Methylphenidate/MTHP	Methylphenidate	300	
Opiates / Morphine/MOP/OPI	Morphine	300	
Opiates 2000/OPI	Morphine	50 / 100 / 200 / 2000	40 / 50 / 300 / 2000
Oxycodone/OXY	Oxycodone	100 / 300	20 / 40 / 50
Phencyclidine/PCP	Phencyclidine	25	10
Propoxyphene/PPX	Propoxyphene	100 / 300	50
Pregabalin/PGB	Pregabalin	500	
Marijuana/THC	Delta-9-THC-COOH	5 / 10 / 15 / 20 / 25 / 40 / 50 / 100 / 200 / 300 / 600	12 / 15 / 20 / 25 / 50
Tramadol/TRA	Tramadol	100 / 200 / 300	50
Tricyclic Antidepressants/TCAs	Nortriptyline	100 / 1000	
Zolpidem/ZOLP	Zolpidem	10	
Zopiclone/ZOP	Zopiclone	50	
6-Monoacetylmorphine/6-MAM	6-Monoacetylmorphine	10	
Carisoprodol/SOMA	Carisoprodol	1000	
Xylazine/XYL	Xylazine	1000	
Alcohol	ALC	0.04%	>0.02%

More items for urine and saliva tests are not included in the table.
The cut off of each item can be customized



“Committed to Health, Driven by Quality”



 www.oryca.com

 +971 4 284 2810

 contact@oryca.com

 Sheikh Saeed Al Maktoum Complex Jebel Ali Industrial Area First, Unit no 198,
Dubai, UAE