

**Prometheus™**  
「Your Trusted Rapid Tests」

# Prometheus

2024 PRODUCT CATALOG



**DRUG OF ABUSE RAPID TEST PAGE 2**

**INFECTIOUS DISEASE TEST CASSETTES PAGE 16**

**VETERINARY TEST CASSETTES PAGE 18**

# THE DRUG TEST CUP LINE

## 510K-APPROVED AND CLIA-WAIVED CUPS FOR DOA URINE TESTS

Our drug test cup product line includes multipanel urine drug test cups designed to be compatible with 30+ CLIA waived & FUD drug options. u-Cups™ are CLIA waived and FDA approved, with many 510K cleared drug panels available. All of these drugs are set to SAMHSA recommended cut-off levels for CLIA waived/510K drugs.

An optional, built-in test is available for up to six different adulterants, including nitrite (NIT), pH, creatinine (CRE), glutaraldehyde (GLUT), specific gravity (S.G.), and oxidant (OX).



### 510K (FDA-Approved) Only

ABBREV.	TEST NAME	CALIBRATOR	CUTOFF
AMP300	Amphetamine	D-Amphetamine	300 ng/mL
COC150	Cocaine	Benzoylcocaine	150 ng/mL
mAMP500	Metamphetamine	D-Metamphetamine	500 ng/mL

### Forensic Use Only (FUD)

ABBREV.	TEST NAME	CALIBRATOR	CUTOFF
ALC	Alcohol	Alcohol	0.04%
APAP	Acetaminophen	Acetaminophen	2000 ng/mL
COT	Cotinine	Cotinine	200 ng/mL
ETG	Ethyl Glucuronide	Ethyl Glucuronide	500 ng/mL
FEN	Fentanyl	Fentanyl	1 ng/mL
K2	K2/Spice	JWH-018/JWH-073	50 ng/mL
KET	Ketamine	Ketamine	1000 ng/mL
LSD	Lysergic Acid Diethylamide	Lysergic Acid Diethylamide	10 ng/mL
MDPV	MDPV	MDPV	300 ng/mL
MEP	Mephedrone	Mephedrone	500 ng/mL
MQL	Methaqualone	Methaqualone	300 ng/mL
RIT	Ritalin	Ritalin	300 ng/mL
SOMA	Carisoprodol	Carisoprodol	1000 ng/mL
TRA	Tramadol	Tramadol	300 ng/mL
ZOLP	Zolpidem	Zolpidem	300 ng/mL
KRA	Kratom	Mitragynine	300 ng/mL 250 ng/mL

### 510K (FDA-Approved) & CLIA-Waived

ABBREV.	TEST NAME	CALIBRATOR	CUTOFF
AMP	Amphetamine	D-Amphetamine	1000 ng/mL
BAR	Barbiturates	Secobarbital	300 ng/mL
BZO	Benzodiazepines	Oxazepam	300 ng/mL
BUP	Buprenorphine	Buprenorphine	10 ng/mL
COC	Cocaine	Benzoylcocaine	300 ng/mL
THC	Marijuana	Delta-9-THC-COOH	50 ng/mL
MTD	Methadone	Methadone	300 ng/mL
mAMP	Metamphetamine	D-Metamphetamine	1000 ng/mL
MDMA	MDMA	MDMA	500 ng/mL
OPI300	Opiates 300	Morphine	300 ng/mL
OPI2000	Opiates 2000	Morphine	2000 ng/mL
OXY	Oxycodone	Oxycodone	100 ng/mL
PCP	Phencyclidine	Phencyclidine	25 ng/mL
PPX	Propoxyphene	Propoxyphene	300 ng/mL
TCA	Tricyclic Antidepressants	Nortriptyline	1000 ng/mL

# ROUND CUP TEST

## AN ADVANCED URINE DRUG SCREENING DEVICE



1. Easy to use, read results within 3-5 minutes with long-lasting color codes.
2. Leak proof with screw-lid locking mechanism.
3. High accuracy
4. Set to SAMHSA cut-off levels.
5. Certifications: CLIA, ISO13485, FDA, 510K

### Round Cup Configurations

PRODUCT ID	DESCRIPTION	DOA TESTS INCLUDED	CERTIFICATIONS
DOW-156-011	Round Cup 5-Panel	AMP, COC, mAMP, THC, OPI2000	CLIA/FDA, CE, ISO13485, FDA Cleared
DOW-166-011	Round Cup 6-Panel	COC, BZO, mAMP, OPI, OXY, THC	CLIA/FDA, CE, ISO13485, FDA Cleared
DOW-1106-011	Round Cup 10-Panel	AMP, BUP, BZO, COC, MDMA, MTD, OPI300, OXY, THC, mAMP	CLIA/FDA, CE, ISO13485, FDA Cleared
DOE-1126-021	Round Cup 12-Panel with Adulteration	AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI300, OXY, PCP, THC + pH, CRE, SG	CLIA/FDA, CE, ISO13485, FDA Cleared
DOW-1126-011	Round Cup 12-Panel	AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI300, OXY, PCP, THC	CLIA/FDA, CE, ISO13485, FDA Cleared
DOW-A1136-011	Round Cup 13-Panel	AMP300, BZO300, COC150, TRA300, ETG500, FEN100, K2-50, MDMA500, MTD300, MOP300, OXY100, THC50, mAMP1000	Forensic use only (FUD)
DOW-A1136-021	Round Cup 13-Panel	AMP300, BZO300, COC150, BUP10, ETG500, FEN100, K2-50, MDMA500, MTD300, MOP300, OXY100, THC50, mAMP1000	Forensic use only (FUD)
DOW-1166-011	Round Cup 16-Panel	FEN, TRA, EtG, OPI, OXY, THC, COC, BZO, AMP, mAMP, BUP, K2, MTD, PCP, MDMA, BAR	Forensic use only (FUD)
DOW-1166-021	Round Cup 16-Panel	FEN, TRA, EtG, MOP, OXY, THC, COC, BZO, AMP, mAMP, BUP, K2, MTD, PCP, MDMA, BAR	Forensic use only (FUD)
U-DOA-3124	Round Cup 13-Panel	MET, COC, MDMA, MOP30, MTD, BZO, BAR, KET, BUP, THC, TRA, K3, PGB, + CRE/S.G./pH	Forensic use only (FUD)
DOX-116-1226	Round Cup CUSTOM Configuration	Select 1 to 22 drug tests plus up to 6 adulterant tests	Depends upon tests selected; see page 2

### ADULTERANT TESTS

Nitrite (NIT) ▪ pH ▪ Creatinine (CRE) ▪ Glutaraldehyde (GLUT) ▪ Specific Gravity (SG) ▪ Oxidant/Pyridinium Chlorochromate (PCC)





# MINI CUP TEST

## SMALL, ECONOMICAL ACCURACY

1. Easy to use, read results within 3-5 minutes with long-lasting color codes.
2. Leak proof with screw-lid locking mechanism.
3. High accuracy
4. Set to SAMHSA cut-off levels.
5. Certifications: CLIA, ISO13485, FDA, 510K



PRODUCT ID	DESCRIPTION	DDA TESTS INCLUDED	CERTIFICATIONS
DFX-1X6-CUSTOM	<b>DOA Mini Cup CUSTOM</b>	Select 1 to 18 drug tests plus up to 6 adulterant tests	Depends upon tests selected; see page 2

# FLAT CUP TEST

## EASILY SCAN, FAX, OR PHOTOCOPY RESULTS

1. Easy to use, read results within 3-5 minutes with long-lasting color codes.
2. Leak proof with screw-lid locking mechanism.
3. High accuracy
4. Set to SAMHSA cut-off levels.
5. Certifications: CLIA, ISO13485, FDA, 510K



PRODUCT ID	DESCRIPTION	DDA TESTS INCLUDED	CERTIFICATIONS
DDX-1X6-CUSTOM	<b>DOA Flat Cup CUSTOM</b>	Select 1 to 18 drug tests plus up to 6 adulterant tests	Depends upon tests selected; see page 2

# KEY CUP TEST

## BETTER SECURE YOUR SAMPLES

1. Easy to use, read results within 3-5 minutes with long-lasting color codes.
2. Leak proof with screw-lid locking mechanism.
3. High accuracy
4. Set to SAMHSA cut-off levels.
5. Certifications: CLIA, ISO13485, FDA, 510K



PRODUCT ID	DESCRIPTION	DDA TESTS INCLUDED	CERTIFICATIONS
DEX-1X6-CUSTOM	<b>DOA Key Cup CUSTOM</b>	Select 1 to 18 drug tests plus up to 6 adulterant tests	Depends upon tests selected; see page 2

# CONE CUP TEST

## DESIGNED FOR FEMALE SAMPLE COLLECTION

1. Easy to use, read results within 3-5 minutes with long-lasting color codes.
2. Leak proof with screw-lid locking mechanism.
3. High accuracy
4. Set to SAMHSA cut-off levels.
5. Certifications: CLIA, ISO13485, FDA, 510K



PRODUCT ID	DESCRIPTION	DDA TESTS INCLUDED	CERTIFICATIONS
DGX-1X6-CUSTOM	<b>DOA Cone Cup CUSTOM</b>	Select 1 to 18 drug tests plus up to 6 adulterant tests	Depends upon tests selected; see page 2



# OEM SERVICES

Let our science support your business.  
Custom labels and configurations available.





# TEST CUP ACCESSORIES


STREAMLINE YOUR SCREENING PROCESS WITH VACUTAINER™



## Vacutainers & Collection Cups:


Vacutainer lids are an optional modification for the further customization of your Cup order. These lids seamlessly transfer portions of a sample into the corresponding Vacutainer tube, providing additional flexibility in sample collection and reading of results.

PRODUCT ID	DESCRIPTION
U-VAC-LID	A sample transfer lid for Vacutainer
U-VAC	A sample transfer vacuum sealed tube



# TEST CARD SCREENING


LEAK-PROOF, HIGH-ACCURACY URINE SCREENING DEVICES



## DESIGN FEATURES:

1. One-step urine drug screening device.
2. Easy to use: simply dip and read.
3. Results within 3-5 minutes with long-lasting color codes.
4. Leak-proof with screw-lid locking mechanism.
5. High-accuracy.
6. Set to SAMHSA cut-off levels.
7. Certifications: CLIA, ISO13485, FDA, 510K
8. Test for multiple drugs simultaneously.
9. Cost effective with minimal waste.
10. Packaging: individual aluminum pouch; 25 per box, 1000 per carton.

Format/pattern/shape: Card-Shaped  
 Strip size: 2.4 mm  
 Specimen: Urine  
 Time to results: 3-5 minutes  
 Service: OEM services available  
 Shelf life/validity: 24 months  
 Certifications: CLIA/510K, CE, ISO-13485, FDA Cleared  
 Storage: between 2-30°C  
 Accuracy: over 99%  
 Manufacture: Prometheus Bio Inc.



Easy to read results:



# MULTI-PANEL TEST CARD

## READY-TO-SHIP TEST CARD CONFIGURATIONS

### Multipanel Test Card Configurations: Ready to Ship

PRODUCT ID	DESCRIPTION	DOA TESTS INCLUDED	CERTIFICATIONS
DOA-3134-081	Drug test card 3-Panel	AMP, COC, mAMP	CLIA/510K, CE, ISO13485, FDA Cleared
DOA-3154-031	Drug Test Card 5-Panel	AMP, COC, mAMP, THC, OPI2000	CLIA/510K, CE, ISO13485, FDA Cleared
DOA-1104-011	Drug Test Card 10-Panel	AMP, BAR, BZO, COC, MTD, OPI2000, OXY, PCP, THC, mAMP	CLIA/510K, CE, ISO13485, FDA Cleared
DOA-1124-241	Drug Test Card 12-Panel	AMP, BAR, BUP, BZO, COC, MDMA, MTD, OPI2000, OXY, PCP, THC, mAMP	CLIA/510K, CE, ISO13485, FDA Cleared
DUE-1124-401	Drug Test Card 12-Panel with Adulteration Strips	AMP, BAR, BUP, BZO, COC, MDMA, MTD, OPI2000, OXY, PCP, THC, mAMP + pH, CRE, SG	CLIA/510K, CE, ISO13485, FDA Cleared

### Single-Panel Test Card Configurations

PRODUCT ID	DESCRIPTION	DOA TEST	CERTIFICATIONS
BAL-114	ALC Test Card	Alcohol (ALC)	Forensic Use Only (FUO)
DAM-114	1-Panel AMP Test Card	Amphetamine (AMP, 1000 ng/mL)	CLIA-waived, FDA cleared
DBA-114	1-Panel BAR Test Card	Barbiturates (BAR, 300 ng/mL)	CLIA-waived, FDA cleared
DBU-114	1-Panel BUP Test Card	Buprenorphine (BUP, 10 ng/mL)	CLIA-waived, FDA cleared
DBZ-114	1-Panel BZO Test Card	Benzodiazepines (BZO, 300 ng/mL)	CLIA-waived, FDA cleared
DCO-114	1-Panel COC Test Card	Cocaine (COC, 300 ng/mL)	CLIA-waived, FDA cleared
DET-114	1-Panel ETG Test Card	Ethyl Glucuronide (ETG, 500 ng/mL)	Forensic Use Only (FUO)
DFE-114	1-Panel FEN Test Card	Fentanyl (FEN, 100 ng/mL)	Forensic Use Only (FUO)
DK2-114	1-Panel K2 Test Card	K2/Spice (K2, 50 ng/mL)	Forensic Use Only (FUO)
DKE-114	1-Panel KET Test Card	Ketamine (KET, 1000 ng/mL)	Forensic Use Only (FUO)
DMD-114	1-Panel MDMA Test Card	MDMA (MDMA, 500 ng/mL)	CLIA-waived, FDA cleared
DMDP-114	1-Panel MDPV Test Card	MDPV (MDPV, 300 ng/mL)	Forensic Use Only (FUO)
DMA-114	1-Panel MET Test Card	Methamphetamine (MET, 1000 ng/mL)	CLIA-waived, FDA cleared
DMO-114	1-Panel MOP Test Card	Morphine (MOP 300 ng/mL)	CLIA-waived, FDA cleared
DMT-114	1-Panel MTD Test Card	Methaqualone (MTD, 300 ng/mL)	CLIA-waived, FDA cleared
DOP-114	1-Panel OPI Test Card	Opiates (OPI, 2000 ng/mL)	CLIA-waived, FDA cleared
DOX-114	1-Panel OXY Test Card	Oxycodone (OXY, 300 ng/mL)	CLIA-waived, FDA cleared
DPC-114	1-Panel PCP Test Card	Phencyclidine (PCP, 25 ng/mL)	CLIA-waived, FDA cleared
DPP-114	1-Panel PPX Test Card	Propoxyphene (PPX, 300 ng/mL)	CLIA-waived, FDA cleared
DTH-114	1-Panel THC Test Card	Marijuana (THC, 50 ng/mL)	CLIA-waived, FDA cleared





# Customize for Your Demographic

SELECT UP TO NINETEENS DRUG TESTS + ADULTERANTS



# DRUG TESTS STRIP

FURTHER CUSTOMIZE YOUR DRUG SCREENING PROCESS

- Test for one particular drug type
- Simply Dip & Read
- 99% accuracy, cost effective & minimum waste
- Set to SAMHSA cut-off levels
- Packaging: 50/Box, 2000/Carton

## 510K (FDA-Approved) & CLIA-Waived

ABBREV.	TEST NAME	CALIBRATOR	CUTOFF
AMP	Amphetamine	D-Amphetamine	1000 ng/mL
BAR	Barbiturates	Secobarbital	300 ng/mL
BZO	Benzodiazepines	Oxazepam	300 ng/mL
BUP	Buprenorphine	Buprenorphine	10 ng/mL
COC	Cocaine	Benzoyllecgonine	300 ng/mL
THC	Marijuana	Delta-9-THC-COOH	50 ng/mL
MTD	Methadone	Methadone	300 ng/mL
mAMP	Metamphetamine	D-Metamphetamine	1000 ng/mL
MDMA	MDMA	MDMA	500 ng/mL
OPI300	Opiates 300	Morphine	300 ng/mL
OPI2000	Opiates 2000	Morphine	2000 ng/mL
OXY	Oxycodone	Oxycodone	100 ng/mL
PCP	Phencyclidine	Phencyclidine	25 ng/mL
PPX	Propoxyphene	Propoxyphene	300 ng/mL
TCA	Tricyclic Antidepressants	Nortriptyline	1000 ng/mL

## Forensic Use Only (FUO)

ABBREV.	TEST NAME	CALIBRATOR	CUTOFF
ALC	Alcohol	Alcohol	0.04%
APAP	Acetaminophen	Acetaminophen	2000 ng/mL
COT	Cotinine	Cotinine	200 ng/mL
ETG	Ethyl Glucuronide	Ethyl Glucuronide	500 ng/mL
FEN	Fentanyl	Fentanyl	1,100 ng/mL
K2	K2/Spice	JWH-018/JWH-073	50 ng/mL
KET	Ketamine	Ketamine	1000 ng/mL
LSD	Lysergic Acid Diethylamide	Lysergic Acid Diethylamide	10 ng/mL
MDPV	MDPV	MDPV	300 ng/mL
MEP	Mephedrone	Mephedrone	500 ng/mL
MQL	Methaqualone	Methaqualone	300 ng/mL
RIT	Ritalin	Ritalin	300 ng/mL
SOMA	Carisoprodol	Carisoprodol	1000 ng/mL
TRA	Tramadol	Tramadol	300 ng/mL
ZOLP	Zolpidem	Zolpidem	300 ng/mL






















## 510K (FDA-Approved) Only

ABBREV.	TEST NAME	CALIBRATOR	CUTOFF
AMP300	Amphetamine	D-Amphetamine	300 ng/mL
COC150	Cocaine	Benzoyllecgonine	150 ng/mL
mAMP500	Metamphetamine	D-Metamphetamine	500 ng/mL














**CUSTOMIZATION TIP:** if one test on a multi-panel Card is Forensic Use Only (FUO), the entire card will be classified as FUO.

PRODUCT ID	DESCRIPTION	DOA TEST	CERTIFICATIONS
BAL-101	ALC Test Strip	Alcohol (ALC)	Forensic Use Only (FUO)
DAM-101	AMP Test Strip	Amphetamine (1000 ng/mL)	FDA Cleared
DAM-A101	AMP300 Test Strip	Amphetamine (300 ng/mL)	FDA Cleared
DAP-A101	APAP5000 Test Strip	Acetaminophen (5000 ng/mL)	Forensic Use Only (FUO)
DBA-101	BAR Test Strip	Barbiturates (300 ng/mL)	Forensic Use Only (FUO)
DBU-101	BUP Test Strip	Buprenorphine (10 ng/mL)	FDA Cleared
DBZ-101	BZO Test Strip	Benzodiazepines (300 ng/mL)	FDA Cleared
DCO-101	COC Test Strip	Cocaine (300 ng/mL)	FDA Cleared
DET-101	ETG Test Strip	Ethyl Glucuronide (500 ng/mL)	Forensic Use Only (FUO)
DFE-101	FEN Test Strip	Fentanyl (100 ng/mL)	Forensic Use Only (FUO)
DK2-50	K2 Test Strip	K2/Spice (50 ng/mL)	Forensic Use Only (FUO)
DKE-101	KET Test Strip	Ketamine (1000 ng/mL)	Forensic Use Only (FUO)
DLS-101	LSD Test Strip	LSD (10 ng/mL)	Forensic Use Only (FUO)
DMD-101	MDMA Test Strip	MDMA (500 ng/mL)	FDA Cleared
DMDP-101	MDPV Test Strip	MDPV (300 ng/mL)	Forensic Use Only (FUO)
U-MQL-114L	MQL Test Strip	Methaqualone (300 ng/mL)	Forensic Use Only (FUO)

# TEST STRIPS CONTINUED

PRODUCT ID	DESCRIPTION	DOA TEST	CERTIFICATIONS	
DCT-101	<b>COT Test Strip</b>	Cotinine (200 ng/mL)	Forensic Use Only (FUO)	
DET-A101	<b>ETG300 Test Strip</b>	ETG (300 ng/mL)	Forensic Use Only (FUO)	
DET-101	<b>ETG Test Strip</b>	ETG (500 ng/mL)	Forensic Use Only (FUO)	
DFE-B101 DEF-101	<b>FEN1 Test Strip</b> <b>FEN Test Strip</b>	Fentanyl (1 ng/mL) Fentanyl (100 ng/mL)	Forensic Use Only (FUO) Forensic Use Only (FUO)	
DK2-101	<b>K2 Test Strip</b>	K2/Spice (50 ng/mL)	Forensic Use Only (FUO)	
DK2-B101	<b>K2-20 Test Strip</b>	K2/Spice (20 ng/mL)	Forensic Use Only (FUO)	
DK2-A101	<b>K2-25 Test Strip</b>	K2/Spice (25 ng/mL)	Forensic Use Only (FUO)	
DKE-101	<b>KET Test Strip</b>	Ketamine (1000 ng/mL)	Forensic Use Only (FUO)	
DLS-101	<b>LSD Test Strip</b>	Lysergic Acid Diethylamide (10 ng/mL)	Forensic Use Only (FUO)	
LSD-A101	<b>LSD3 Test Strip</b>	Lysergic Acid Diethylamide (3 ng/mL)	Forensic Use Only (FUO)	
DMA-101	<b>mAMP Test Strip</b>	Methamphetamine (1000 ng/mL)	FDA Cleared	
DMA-C101 DMA-B101	<b>mAMP200 Test Strip</b> <b>mAMP300 Test Strip</b>	Methamphetamine (200 ng/mL) Methamphetamine (300 ng/mL)	Forensic Use Only (FUO) Forensic Use Only (FUO)	
DMA-D101	<b>mAMP50 Test Strip</b>	Methamphetamine (50 ng/mL)	Forensic Use Only (FUO)	
DMA-A101	<b>mAMP500 Test Strip</b>	Methamphetamine (500 ng/mL)	Forensic Use Only (FUO)	
DMD-101	<b>MDMA Test Strip</b>	MDMA (500 ng/mL)	FDA Cleared	
DMDP-101	<b>MDPV Test Strip</b>	Methylenedioxypropylvalerone (300 ng/mL)	Forensic Use Only (FUO)	
DMDP-A101	<b>MDPV500 Test Strip</b>	Methylenedioxypropylvalerone (500 ng/mL)	Forensic Use Only (FUO)	
DME-101	<b>MEP Test Strip</b>	Mephedrone (500 ng/mL)	Forensic Use Only (FUO)	
DMQ-101	<b>MQL Test Strip</b>	Methaqualone (300 ng/mL)	Forensic Use Only (FUO)	
DMT-101	<b>MTD Test Strip</b>	Methadone (300 ng/mL)	FDA Cleared	
DOP-B101	<b>OPI100 Test Strip</b>	Opiates (100 ng/mL)	Forensic Use Only (FUO)	

# TEST STRIPS CONTINUED

PRODUCT ID	DESCRIPTION	DOA TEST	CERTIFICATIONS	
DOP-101	<b>OPI Test Strip</b>	Opiates (2000 ng/mL)	FDA Cleared	
DOP-D101	<b>OPI300 Test Strip</b>	Opiates (300 ng/mL)	FDA Cleared	
DOP-C101	<b>OPI150 Test Strip</b>	Opiates (150 ng/mL)	Forensic Use Only (FUO)	
DOX-101	<b>OXY Test Strip</b>	Oxycodone (100 ng/mL)	FDA Cleared	
DOX-A101	<b>OXY300 Test Strip</b>	Oxycodone (300 ng/mL)	Forensic Use Only (FUO)	
DPC-101	<b>PCP Test Strip</b>	Phencyclidine (25 ng/mL)	FDA Cleared	
DPP-101	<b>PPX Test Strip</b>	Propoxyphene (300 ng/mL)	FDA Cleared	
DRI-101	<b>RIT Test Strip</b>	Ritalin (300 ng/mL)	Forensic Use Only (FUO)	
DSO-101	<b>SOMA Test Strip</b>	Carisoprodol (1000 ng/mL)	Forensic Use Only (FUO)	
DTC-101	<b>TCA Test Strip</b>	Tricyclic Antidepressants (1000 ng/mL)	FDA Cleared	
DTH-101	<b>THC Test Strip</b>	Marijuana (50 ng/mL)	FDA Cleared	
DTH-H101	<b>THC10 Test Strip</b>	Marijuana (10 ng/mL)	Forensic Use Only (FUO)	
DTH-K101	<b>THC100 Test Strip</b>	Marijuana (100 ng/mL)	Forensic Use Only (FUO)	
DTH-B101	<b>THC20 Test Strip</b>	Marijuana (20 ng/mL)	Forensic Use Only (FUO)	
DTH-J101	<b>THC200 Test Strip</b>	Marijuana (200 ng/mL)	Forensic Use Only (FUO)	
DTH-G101	<b>THC25 Test Strip</b>	Marijuana (25 ng/mL)	Forensic Use Only (FUO)	
DTH-A101	<b>THC40 Test Strip</b>	Marijuana (40 ng/mL)	Forensic Use Only (FUO)	
DTH-I101	<b>THC50 Test Strip</b>	Marijuana (50 ng/mL)	Forensic Use Only (FUO)	
DTR-A101	<b>TRA100 Test Strip</b>	Tramadol (100 ng/mL)	Forensic Use Only (FUO)	
DTR-B101 DTR-101	<b>TRA200 Test Strip</b> <b>TRA Test Strip</b>	Tramadol (200 ng/mL) Tramadol (300 ng/mL)	Forensic Use Only (FUO) Forensic Use Only (FUO)	
DZO-101	<b>ZOLP Test Strip</b>	Zolpidem (10 ng/mL)	Forensic Use Only (FUO)	



# SALIVA ROUND CUP

## ECONOMICAL AND LESS INVASIVE THAN URINE

1. Easy to use, read results within 3-5 minutes with long lasting color codes.
2. Small size significantly saves on shipping cost
3. Non-invasive and decent collection
4. Three step procedure: 1) Swab 2) Collect 3) Read
5. High-accuracy
6. More sensitive than urine drug testing
7. Minimal chances of adulteration and contamination
8. Avoids cross gender observation
9. Set to SAMHSA cut-off levels
10. Forensic use only (FUO)



PRODUCT ID	DESCRIPTION	DOA TESTS INCLUDED	CERTIFICATIONS
DSC-CUSTOM	SALIVA ROUND CUP CUSTOM	Select 1 to 12 drug tests	Forensic use only (FUO)

# SALIVA TEST CUBE

## ECONOMICAL AND LESS INVASIVE THAN URINE

1. Easy to use, read results within 3-5 minutes with long lasting color codes.
2. Small size significantly saves on shipping cost
3. Non-invasive and decent collection
4. Three step procedure: 1) Swab 2) Collect 3) Read
5. High-accuracy
6. More sensitive than urine drug testing
7. Minimal chances of adulteration and contamination
8. Avoids cross gender observation
9. Set to SAMHSA cut-off levels
10. Forensic use only (FUO)



PRODUCT ID	DESCRIPTION	DOA TESTS INCLUDED	CERTIFICATIONS
DSE-254	Saliva Cube 5-Panel	AMP, COC, mAMP, THC, OPI2000	Forensic use only (FUO)
DSE-7126-011	Saliva Cube 12-Panel	ALC, AMP, BAR, BUP, BZO, COC, MTD, OPI2000, OXY, PCP, THC, mAMP	Forensic use only (FUO)
DSE-CUSTOM	Saliva Cube CUSTOM	Select 1 to 14 drug tests	Forensic use only (FUO)

# CUP READER


## MORE ACCURATE THAN THE HUMAN EYE...AND FASTER, TOO.

1. Compact, standalone device
2. Reads round UCP cups
3. Accurate, reproducible readout
4. Minimum User Operation
5. Quick Scan
6. Optional PC link for data export



## Why use the cup Reader?

Our Cup Reader features an intuitive user interface and powerful processing capabilities. It automatically recognizes which UCP cup has been inserted and how many Drug of Abuse (DOA) screening strips it has, and will complete scanning all of them within 1 minute.



SAVE TO USB WITH ENCRYPTION OPTIONS OR DIRECTLY TO COMPUTER WITH OPTIONAL LINK.

# SALIVA TEST MIDSTREAM

## MORE CONCISE AND CONVENIENT

1. Easy to use, just hold the saliva sponge tip in your mouth and wait for the colored lines to appear to get your results
2. Minimal chances of adulteration and contamination
3. Non-invasive and decent collection
4. Two step procedure: 1) Sample collect 2) Read result
5. High-accuracy





# INFECTIOUS DISEASE LINE

## A METHOD TO ASSIST IN THE DIAGNOSIS OF INFECTIOUS DISEASES

Our infectious disease product line includes Respiratory Pathogens Test Series, Vector Borne Pathogens Test Series, Digestive System Pathogens Test Series, Maternal and child health Test Series and other infectious disease test. All of these products are CE approved.

These test can auxiliary diagnosis, treatment, prevention, outbreak investigation, surveillance, and vaccine development. Its timely and accurate implementation contributes to better patient outcomes, disease control, and public health strategies.



### Vector Borne Pathogens Test

ABBREV.	TEST NAME	SPECIMEN	TEST/KIT
DV IgG/IgM	Dengue IgG/IgM	Whole blood/Serum /Plasma	1, 10, 25, 40
DV NS1	Dengue NS1	Whole blood/Serum /Plasma	1, 10, 25, 40
DV NS1+ IgG/IgM	Dengue NS1+IgG/IgM	Whole blood/Serum /Plasma	1, 10, 25, 40
Pf/Pv	Malaria Pf/Pv	Whole blood	1, 10, 25, 40
Pf/Pan	Malaria Pf/Pan	Whole blood	1, 10, 25, 40
Pf	Malaria Pf	Whole blood	1, 10, 25, 40

### Maternal and child health Test

ABBREV.	TEST NAME	SPECIMEN	TEST/KIT
CT-IgM	Chlamydia trachomatis IgM	Whole blood/Serum /Plasma	1, 10, 25, 40
CT-Ag	Chlamydia trachomatis Ag	Secretion	1, 10, 25, 40
HSV-IgG/IgM	Herpes simplex Virus IgG/IgM	Whole blood/Serum /Plasma	1, 10, 25, 40
GBS-Ag	Group B Streptococcus Antigen	Secretion	1, 10, 25, 40
TOXO-IgM	Toxoplasma gondii IgM Antibody	Whole blood/Serum /Plasma	1, 10, 25, 40
TOXO-IgG	Toxoplasma gondii IgG Antibody	Whole blood/Serum /Plasma	1, 10, 25, 40

### Respiratory Pathogens Test

ABBREV.	TEST NAME	SPECIMEN	TEST/KIT
Flu A/B	Influenza A/B Ag	Secretion	1, 10, 25, 40
4 in 1	SARS-CoV-2/Flu A/B /RSV-Ag 4in1	Secretion	1, 10, 25, 40
3 in 1	SARS-CoV-2/Flu A/B -Ag 3in1	Secretion	1, 10, 25, 40
MP-IgM	Mycoplasma pneumoniae IgM	Whole blood/Serum /Plasma	1, 10, 25, 40
Cpn-IgM	Chlamydia pneumoniae IgM	Whole blood/Serum /Plasma	1, 10, 25, 40
TB-IgG/IgM	Mycobacterium tuberculosis IgG/IgM	Whole blood/Serum /Plasma	1, 10, 25, 40
TB-Ab	Mycobacterium tuberculosis Ab	Whole blood/Serum /Plasma	1, 10, 25, 40
GAS-Ag	Group A Streptococcus Ag	Secretion	1, 10, 25, 40
RSV-Ag	Respiratory syncytial Virus Ag	Secretion	1, 10, 25, 40
SARS-COV-2 Ag	SARS-COV-2 Antigen	Secretion	1, 10, 25, 40
COVID-19 IgM/IgG	COVID-19 IgM/IgG	Whole blood/Serum /Plasma	1, 10, 25, 40
COVID-19 NAb	COVID-19 Neutralization Antibody	Whole blood/Serum /Plasma	1, 10, 25, 40

# DIGESTIVE SYSTEM PATHOGENS

## IDENTIFICATION AND DETECTION OF MICROORGANISMS

1. Quick results: Results are available in just 5-20 minutes (Depends on the product), making them ideal for use in various settings.
2. Simple and convenient: micro sample detection.
3. Transport Easy: the kits are stored at room temperature and do not require refrigeration or cold chain transport.
4. Shelf life: 12-24 month long shelf life and can be stored easily.
5. Complete qualifications: CE, ISO13485 and NMPA.




Digestive System Pathogens Test Devices are used to detect the related pathogens of gastrointestinal infection; assist in the diagnosis and treatment of foodborne diseases; distinguish digestive system bacteria and virus infections.

### DSP Test Configurations

ABBREV.	TEST NAME	SPECIMEN	TEST/KIT
HAV-IgM	Hepatitis A virus IgM Antibodies	Whole blood/Serum/Plasma	1, 10, 25, 40
HP-Ab	Helicobacter pylori Antibodies	Whole blood/Serum/Plasma	1, 10, 25, 40
HP-Ag	Helicobacter pylori Antigen	Feces	1, 10, 25, 40
Sa-IgG/IgM	Salmonella typhi IgG/IgM Antibodies	Whole blood/Serum/Plasma	1, 10, 25, 40
Sa-Ag	Salmonella typhi Antigen	Feces	1, 10, 25, 40
AdV-Ag	Adenovirus Antigen	Feces	1, 10, 25, 40
RoV-Ag	Rotavirus Antigen	Feces	1, 10, 25, 40
EV71-IgM	Enterovirus 71 IgM Antibodies	Whole blood/Serum/Plasma	1, 10, 25, 40
COXA16-IgM	Coxsackie virus type A16 IgM Antibody	Whole blood/Serum/Plasma	1, 10, 25, 40

### OTHERS ID TESTS

Human Brucella IgM ▪ Clostridium difficile Toxin A/B Ag  
Fecal Occult Blood ▪ Monkeypox Virus IgM/IgG  
Monkeypox Virus Antigen



# CANINE RAPID TESTS


## TEST KITS FOR USE WITHOUT ANY INSTRUMENT

These kits are used for the in vitro qualitatively detection of disease in canine secretions, feces or vomit samples.

Canine rapid tests are medical devices used to quickly diagnose and monitor dogs health conditions. Included heartworm infection, canine parvovirus, distemper, Lyme disease, etc. The tests are easy to perform and provide rapid and accurate results. Canine rapid tests are commonly used in veterinary clinics, animal shelters, and rescue organizations.

Early detection of these diseases is essential for timely treatment and preventing the spread of infection/treatment, prevention, outbreak investigation.

ABBREV.	TEST NAME	SPECIMEN	TEST/KIT
CPV Ag	Canine Parvovirus Antigen	Feces/Vomit	1, 10
CCV Ag	Canine Coronavirus Antigen	Feces/Vomit	1, 10
CPV + CCV Ag Combo	CPV + CCV Ag Combo	Feces/Vomit	1, 10
CIV Ag	Canine Influenza Virus Antigen	Secretions	1, 10
CDV Ag	Canine Distemper Virus Antigen	Secretions	1, 10
CAV Ag	Canine Adenovirus Type - II Antigen	Secretions	1, 10
CDV+CAV+CIV-Ag	CDV+CAV-II+CIV Antigen Combo	Secretions	1, 10
TOXO Ag	Canine Toxoplasma Antigen	Whole blood/Serum/Plasma	1, 10
<b>Storage</b>	Cool, dark and dry place at 2-30°C (Do not freeze)		
<b>Shelf Life/Validity</b>	18 Months		
<b>Time to Result</b>	The results will be get within 5-10 minutes		
<b>Packaging Specifications</b>	1 kits/box, 10 kits/box		



# FELINE RAPID TESTS

## TEST KITS FOR USE WITHOUT ANY INSTRUMENT

These kits are used for the in vitro qualitatively detection of disease in feline feces, vomit, ascites, hydrothorax, eye secretion or secretions samples.

Feline rapid tests are medical devices used to quickly diagnose and monitor cats health conditions. Included feline panleucopenia virus (FPV), feline calicivirus (FCV), feline herpes virus(FHV), etc. The tests are easy to perform and provide rapid and accurate results. Feline rapid tests are commonly used in veterinary clinics, animals helters, and rescue organizations.

Early detection of these diseases is essential for timely treatment and preventing the spread of infection/treatment, prevention, outbreak investigation.

ABBREV.	TEST NAME	SPECIMEN	TEST/KIT
FPV Ag	Feline Herpes Virus Antigen	Feces/Vomit	1, 10
FCoV Ag	Feline Coronavirus Antigen	Ascites/Hydrothorax	1, 10
FPV + FCoV Ag Combo	FPV + FCoV Ag Combo	Feces	1, 10
FHV Ag	Feline Herpes Virus Antigen	Eye Secretion	1, 10
FCV Ag	Feline Calicivirus Antigen	Secretions	1, 10
TOXO Ag	Feline Toxoplasma Antigen	Feces	1, 10
<b>Storage</b>	Cool, dark and dry place at 2-30°C (Do not freeze)		
<b>Shelf Life/Validity</b>	18 Months		
<b>Time to Result</b>	The results will be get within 5-10 minutes		
<b>Packaging Specifications</b>	1 kits/box, 10 kits/box		
<b>Strip Width Size</b>	The regular strip width is 3mm or 4 mm		
<b>Principle</b>	Colloidal gold method		







Sheikh Saeed Al Maktoum Complex, Jebel Ali Industrial Area First, Unit no: 198 Dubai, UAE  
+971(4)2842810                      [contact@oryca.com](mailto:contact@oryca.com)                      [www.oryca.com](http://www.oryca.com)